


Red Light Therapy Panel

Glow Light Series - GL960

This vertical red light therapy panel features touchscreen control and seven combined wavelengths, 480, 630, 660, 810, 830, 850, and 940 nm, providing hands-free support for skin health, muscle recovery, joint comfort, and overall wellness.

Key Data at a Glance

Power	LED Quantity
1456W	960 pcs
Wavelengths	Irradiance at 6 inches
7 Wavelengths	190mw/cm ²
Color	Feature
Black and White	OEM/ODM available



Use Note

Start with shorter sessions and increase gradually based on individual tolerance. Avoid looking directly into the LEDs. Use eye protection for comfort, especially during face or upper-body sessions.

Product Specifications

Glow Light Series - GL960 / Red Light Therapy Panel

Model GL960	Power 1456W	Lens Beam Angle 30 Degree
Irradiance at 6 inches 190mw/cm²	Lifespan >100000 Hours	LED Quantity 960 pcs
Input Voltage 100v-250v	Material Steel	N.W.(KG) 28 Kg
Flicker No Flicker	G.W.(KG) 35 Kg	Unit Package Dimension 195*58*15cm
Working Modes 1. Red/Nir/Blue 2. Any combination	EMF Emission 0µT @ 3inches	Sound(dB)@15cm 52 DB
Cooling Fans Quantity 12pcs	Working Temperature -20°C~40°C	Warranty 2 Years
Product size 1800*473*60mm (70.87*18.62*2.36")	Timer 1-30Min	Remote Control Yes
Wavelength (480/630/660/810/830/850/940nm) customized wavelength		



Panel Configuration

The GL960 supports touch screen control, seven wavelengths, and flexible working modes including Red/NIR/Blue and custom combinations. Its large vertical format is suitable for broader treatment coverage across the back, legs, and full-body wellness routines.

Package Contents

Standard accessories supplied with the GL960 panel.



Power Cable *1



Mounting Hardware *1



Lifting Rope *1



Goggles *1



Controller *1



User Manual *1



Door Hanging Kit *1

Additional Treatment Tips

● Session Duration & Frequency

For deep tissue treatment, use 15-30 min sessions for each target area at approx 15-30cm from the device, with a max of two to three treatment areas per day. For skin or superficial therapy, a closer distance of around 5-15 cm with 5-10 minutes per session is usually suitable. Time and distance may vary by individual tolerance.

● Eye Protection

Avoid looking directly into the LEDs, as the light can be uncomfortable. Goggles are included in each package for comfort and added protection, especially when treating areas near the face.

● Clothing

Make sure the treatment area is fully exposed to the light and not covered by fabric or clothing. Areas that are not being treated may be covered if preferred.

● Effects

It is common to feel gentle warmth during and after a session. With consistent use, users may notice muscle relaxation, reduced stiffness, and improved comfort. If discomfort or excessive sensitivity occurs, reduce exposure time or increase the distance from the light.

Frequently Asked Questions

Quick answers for customer use, installation, and basic safety review.

● **1. Do you need glasses for red light therapy?**

Yes. Eye protection is recommended, especially with bright panels. Do not stare directly into the LEDs. Use the provided goggles or keep your eyes closed during treatment.

● **2. How to use red light therapy for sleep?**

Use it 30-60 minutes before bed for about 10-20 minutes. Keep a comfortable distance, relax during the session, and avoid staring directly at the light.

● **3. What are benefits of red light therapy?**

Red light therapy may support collagen production, skin repair, inflammation balance, muscle recovery, joint comfort, wound-healing support, hair wellness, and overall cellular energy.

● **4. Who should not use red light therapy?**

People with photosensitive conditions, active cancer, epilepsy, pregnancy, open lesions, or photosensitizing medications should consult a healthcare professional before using red light therapy.

● **5. Does red light therapy help with sunburn?**

Red light therapy may help calm redness, support skin repair, and ease discomfort after sun exposure. It should not replace medical care for severe sunburn or blistering.

● **6. Is red light therapy a hype?**

Red light therapy is real science, but often over-marketed. It may support skin, recovery, inflammation balance, and hair wellness, but it is not a cure-all or instant-result treatment.

● **7. Can red light therapy help with beards and eyebrows too?**

It may support beard and eyebrow growth by improving local circulation, follicle activity, and the growth environment. Evidence is stronger for scalp hair, and results may vary.

● **8. How long should you do red light therapy?**

Most users apply red light therapy for 10-20 minutes per treatment area, 3-5 times per week. Consistency matters more than long sessions, and results usually take several weeks.

● **9. Does red light therapy tan you?**

No. Red light therapy does not use UV rays and does not tan the skin. It works through photobiomodulation, supporting cellular energy, collagen activity, and recovery.

● **10. What is the most important advice for red light therapy users?**

Be patient and consistent. Follow the recommended time, distance, and frequency for your device. More time does not always mean better results; the right dose matters most.

Reference Note

This document is provided for product reference and customer download. It is not a medical document and should not replace professional medical advice.

Shenzhen Medlight Technology Co.,Ltd



Room 306, Bld D, Dingfeng Scientific Park, Dalang Town, Longhua District, Shenzhen, China 518109



Web: www.medledlight.com



Cell: +86-13554796007



Email: sales@medledlight.com