

Red Light Therapy Bed

M7 | M7+

A full-body red and near-infrared light therapy system designed to support skin rejuvenation, muscle recovery, joint comfort and circulation. It has the familiar form of a tanning bed, but uses LED photobiomodulation and emits no ultraviolet light.



OPEN CONFIGURATION



CLOSED CONFIGURATION

HOW IT WORKS & BENEFITS

1 Cellular Energy

Red and near-infrared light is absorbed by mitochondria, supporting ATP production through photobiomodulation.

2 Wavelength Depth

633 and 660 nm target surface tissue; 810, 850 and 940 nm reach deeper muscle and joint structures.

3 Full-Body Coverage

The bed format provides uniform, 360-degree exposure in one automated session.

4 Recovery Support

Designed to support circulation, muscle recovery, joint mobility and inflammation management.

HOW TO USE

FREQUENCY

Do sessions 3 to 5 times per week.

DURATION

Keep each session between 10 and 20 minutes.



PREPARATION

Remove clothing so the light hits your skin directly.

Specifications

PARAMETER	SPECIFICATION
GENERAL INFORMATION	
Model	M7 M7+
Type	Whole Body RED Light Therapy Bed
Body Colors	Lunar White and Crystal Black (custom edition)
LIGHT THERAPY SYSTEM	
Light Source	Taiwan EPISTAR 0.2 W special LED
LED Quantity	13,440 pcs 26,880 pcs
Wavelengths	633 nm, 660 nm, 810 nm, 850 nm, 940 nm
Irradiance	48 to 100 mW/cm ²
Independent Wavelength Control	Yes, independent control of each wavelength
Pulse System	Yes, 0 to 1,000 Hz pulse modulation
LED Lifespan	100,000 hours
CONTROL AND FEATURES	
Control System	Smart Controller 2.0 or Wireless Pad Controller 2.0 (optional)
Heat Abstractor	Self-circulating forced-air cooling system
Wireless PAD Remote	Yes
Bluetooth Stereo Music	Yes
Eyewear Goggles	Yes, included
Power Plug	No
Movable Wheels	Yes
ELECTRICAL SPECIFICATIONS	
Input Voltage	220 V (or 380 V)
Current (220 V)	12 A 21 A
Power (220 V)	2,200 W 4,500 W
DIMENSIONS AND WEIGHT	
Dimensions	2145 mm x 1060 mm x 910 mm (L x W x H)
Inner Height	370 mm
Weight	270 kg
APPLICATIONS	
Usage	Commercial business/spa or home use

Product Dimensions

OVERALL**2145 x 1060 x 910****OPEN HEIGHT****1435****INNER WIDTH****810****INNER HEIGHT****370****END VIEW - OPEN****THREE-QUARTER VIEW - CLOSED****Reference Note**

This document is provided for product reference and customer download. It is not a medical document and should not replace professional medical advice.

Shenzhen Medlight Technology Co., Ltd



Room 306, Bld D, Dingfeng Scientific Park, Dalang Town, Longhua District, Shenzhen, China 518109



Web: www.medledlight.com



Cell: +86-13554796007



Email: sales@medledlight.com